INSTRUCTIONS FOR USE

OMNitech® Titanium Dental Implant System; OMNICON ACTIVE System

EN

1. DESCRIPTION

OMNitech® Titanium Dental Implant System and its sub-model OMNICON ACTIVE products are biocompatible titanium -Gr5 (Ti6Al4V)-made, small size dental implants and prosthetic parts that have been applied to upper and lower jaw bones for the purposes of recovery of missing tooth/teeth, recovery of tooth and jaw functions, and supporting the jaw bones and palate.

Subjected products are single-use and invasive devices. They do not include blood, human or animal tissue, phthalate. Any software or accessory also not included in the products. They are launched in sterile form (gamma sterilization) and do not cause any radioactive or ray scattering. The products are implant devices that directly contacted to the bone/tissue and the duration of the contact is more than 30 days.

OMNitech® Dental Implants: Indicated for single or multiple recovery of teeth and supporting of jaw bones and palate. 3.5 mm platform (diameters 3.7, 4.2, 4.7 mm); and 4.5 mm platform (5.2 mm diameter) implants can be used in the upper and lower posterior. 3.0 platform (diameter 3.2 mm) can be used in the upper and lower anterior.

OMNICON ACTIVE Dental Implants: A model of Omnitech® Dental Implants. Connection of the implants is conical. Indicated for single or multiple recovery of teeth and supporting of jaw bones and palate. 3.4 mm platform (diameters 4.3, 5.0 mm) implants can be used in the upper and lower posterior. 3.0 platform (diameters 3.2, 3.7 mm) can be used in the upper and lower anterior.

Package:

Sterile implant vials contain the following items: an implant body, a transfer part which carries the implant and preventing it to touch inside of the bottle, and a closing screw which is intended to use after screwing the implant. Vials are put into safety boxes individually which once it opens, safety seal junction is teared. Teared seal boxes must not be used.

2. INTENDED USAGE

The main intended use of **OMNitech® Titanium Dental Implant System** is to provide recovery to missing teeth/tooth, missed tooth and jaw functions, and to support the jaw bones and palate by promoting osseointegration of implants with upper and lower jaw bones.

3. INTENDED PERFORMANCE

Intended performance of the products is to provide recovery of missing tooth/teeth, recovery of tooth and jaw functions, and supporting the jaw bones and palate.

4. INDICATIONS

Indications and contra-indications for the system is precisely the same points. The system has the common indications and contra-indications because there are only size differences exists among the products.

- Recovery of missing tooth/teeth.
- Recovery of missed tooth and jaw functions.
- Supporting the jaw bones and palate.

5. COMPATIBILITY

OMNitech® Titanium Dental Implant System surgically compatible with **OMNitech® Surgical Instruments**.

6. CONTRAINDICATIONS

- Vascular conditions
- Uncontrolled diabetes
- Clotting disorders
- Anticoagulant therapy
- Metabolic bone disease
- Chemotherapy or radiation therapy
- Chronic periodontal inflammation
- Insufficient soft tissue coverage
- Metabolic or systemic disorders associated with wound and/or bone healing
- Use of pharmaceuticals that inhibit or alter natural bone remodeling
- Any disorders which inhibit a patient's ability to maintain adequate oral hygiene.
- Foreign body reaction
- Immunodeficiency disorders
- Bleeding disorders

Oral Contraindications:

Include the followings but not limited to:

- Uncontrolled parafunctional habits (e.g. bruxing, clenching, gnawing)
- Insufficient height and/or width of bone
- Insufficient inter-arch space
- Intraoral infection
- Poor or noncompliant patient oral hygiene

7. SURFACES

Implants are made of titanium alloy (Grade 5) (Ti6Al4V). The surface of the implants is created by using TCP particle to create a medium-rough texture. The surfaces are achieved via sandblasting procedure.

8. STERILE PACKAGING

 Implants and included materials (transfer part and closing screw) are marketed **STERILE** by gamma irradiation. Do not re-sterilize.

9. STORAGE & HANDLING

- Do NOT use if packaging is damaged.
- Implants must be stored in a dry place, at room temperature, in their original packaging.
- Dental implants are provided mounted to transfer parts in sterile vials. Transfer part is intended to be used to transport the implant to the prepared surgical site. Appropriate sterilized surgical tools, with retention, are used to transport and place the transfer part mounted implants.
- Do NOT handle implant surfaces directly.
- Users should inspect boxes visually to ensure seals and contents are intact and in their original packaging prior to use.

10. SINGLE-USE PACKAGING

Implants are NOT PERMITTED to be RE-USED.

- The effects of multiple use have not been clinically or scientifically tested, and there is no data to support that the device will perform as designed after multiple uses.
- Reprocessing and resterilization is NOT allowed.
- There is a high risk of infection, unacceptable performance, or device failure if they are not used as specified.

11. SHELF-LIFE

YYYY-MM: Dental implants are considered sterile for five years from the date of initial sterilization. The product expiration date is indicated by the hourglass symbol on the product labels (both inner and outer), followed by year and month of expiration.

12. WARNINGS

- Implants have not been evaluated in terms of safety and compatibility in MRI environment. They have not been tested for heating, migration, image artifact in the MR environment. The safety of implants in the MR environment is unknown.
- Do NOT use if packaging is damaged.
- Implants are NOT PERMITTED to be RE-USED.
- Implants in teared seal boxes must NOT be used.
- Small diameter and angled abutments are not recommended for the posterior region of the mouth.
- **OMNitech**® model 3.0 platform implants shall be used with 3.0 platform prosthetic components. 3.5 platform implants shall be used with 3.5 platform prosthetic components. 4.5 platform implants shall be used with 4.5 platform prosthetic components.
- **OMNICON ACTIVE** sub-model 3.0 platform conical implants shall be used with 3.0 platform prosthetic components. 3.4 platform conical implants shall be used with 3.4 platform prosthetic components.
- Due to conductive nature of metallic implants, electrosurgery is contraindicated.

13. PRECAUTIONS

- For short implants, operators should closely monitor patient for any of the following conditions: peri-implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length.
- If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal.

14. POSSIBLE COMPLICATIONS & ADVERSE REACTIONS

- Loosening, breakage, or falling of implant caused by non-osseointegration, osteoporosis, inhibited revascularization, bone resorption, lack of adhesion, poor bone formation, or surgical failure
- Non-union, malunion, or delayed union of the implant with the bone resulting from a series of causes including smoking, poor bone formation, surgical mistakes etc.
- Pain, discomfort, or abnormal sensitization due to presence of the implant
- Necrosis of bone or tissues
- Infection, nerve damage, sinus damage
- Early or late infection
- Inflammation
- Pyrogenic effects

15. PROSTHETIC COMPONENTS

OMNitech® Titanium Dental Implant System and its sub-model OMNICON ACTIVE prosthetic components are designed to be used in support of a dental implant(s) to provide support for prosthetic restorations in edentulous patients. Prosthetic components are supplied in non-sterile form.

Sterilization of Prosthetic Components: Steam sterilize using a gravity displacement cycle for 30 minutes at 134 °C (270°F), 15-20 psi, and dry time: 15 minutes wrapped in a double pouch.

16. SURGICAL TECHNIQUES FOR IMPLANT PLACEMENT

1. *Pre-Operative Treatment Planning:* During the pre-operative stage, availability of bone-height and width must be determined. Appropriate radiography should be used to determine bone availability, optimal implant location, and to avoid structures such as the mandibular canal, maxillary sinuses, and adjacent teeth.
2. *Surgical Site Preparation:* Follow the corresponding drilling sequence of soft and hard bone preparation respectively.
3. *Insertion Procedure:*
 - a. *Patient Information:* Catalog numbers **REF** are specified on labels. Peel off removable label from the vial and place the label in patient's chart.
 - b. *Opening Outer Box:* Open the cap of the box. The safety seal (extent of outer label) is detached. Take out the vial from the box.
 - c. *Opening the Vial:* Remove vial cap by turning counterclockwise.
 - d. *Remove from Vial:* Mount the appropriately sterilized hex ratchet instrument onto the transfer part-implant assembly by using ratchet. Take the assembly out from the vial with plastic parts and easily remove the plastic parts.
 - e. *Implant to Site:* Carry the implant to the prepared surgical site and thread it for two full rounds. Then remove the transfer part from implant.
 - f. *Implant Insertion:* Continue inserting the implant using appropriate instrumentation. Torque values have great importance in the insertion procedure.
 - g. *Closing Screw:* Use the appropriate instrumentation and unthread the implant closing screw from the plastic part located beneath the implant in the vial and thread it into the implant.
 - h. *Post-Operative Care:* It is recommended that patients use a suitable mouth rinse and perform regular oral hygiene following surgery.
 - i. *Healing Time:* Usually implants are waited to heal for a period of two to four months prior to being restored depending upon bone quality, type, or any compromising medical condition.
 - j. *Instructions:* Users are advised to consult packaged instructions for use for assistance.

SYMBOLS:

Manufacturer

Use-by date

CE marked devices complies with 93/42/EEC requirements.

Sterilized using irradiation

Do not re-use

Do not re-sterilize

Do not use if package is damaged

Single sterile barrier system with protective packaging outside

Consult instructions for use (IFU)

Caution

Catalogue number

Batch code

Medical device

Unique identification number

Non-sterile (for prosthetic components)

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